



Equibiomedic

Venta, alquiler y reparación de equipos biomédicos.

CATÁLOGO

Equipos Biomédicos

Medellín - Colombia

HugeMed Reusable Video Laryngoscope

VL4R



The image quality

- 4inch full-angle LCD high resolution display, Resolution ratio: 800*480RGB
- Rotation: up/down:140° left/right:180°



Functional

- one-click button for photography and video recording
- built-in 8G storage card



Model is complete

- 8 models of blades are designed for various population

HugeMed Disposable Video Laryngoscope

VL4D



Convenient for teaching

- Through the USB data interface, charging, video copy the output
- Extensible external HD interface display



High resolution wide-angle

- 2 megapixels camera
- Field angle ≥ 60 degrees, and the lens angle is between 50° - 60°



Model is complete

- 5 models of blades are designed for various population

Safety

Anti-Fog Lens
Anti-Microbial Handle

Disposable blades
Sterile package
Single use design avoid Cross Infection

Reusable blades

Durable 316 Medical stainless steel
IPX8 waterproof grade design
Support immersion disinfection

Faster

Lightweight, Portable
Faster Installation and Intubation
Faster Removal of blades or Drop blades

High-Definition

The 4 inch anti-reflective display with full view has the high performance visualization capabilities for every intubation
Better View: 2.0 Megapixels Cameras

Easier

Rechargeable li-ion battery, One Shot Button for Recording

Cost Effective


Low cost disposable blades
316 medical stainless reusable blades can be used 1000 times

Technical Specifications

Display	Size: 4inch Full view
	Monitor: LCD
	Resolution Ratio: 800*480RGB
	Aspect Ratio: 4:3
	Video refresh Rata: 30FPS
Data output	Easy to establish and store files
Image output	USB interface
Laryngoscope blade Camera Include	Resolution Ratio 2.0M Pixels
	Field Angle: 66°
	Illuminance ≥800LUX
Handle	Confortable ergonomic short handle design
Battery	Type Rechargeable Lithium Battery
	Long working time ≥120-150 Minutes
	Voltage 3.7V
	Capacity 3200mAh
	Battery life Cycle: More Than 300 times
Charging times: less than 4-6hrs	



ACCESSORIES



MOBILE STAND

A stand-alone configuration for use in a variety of hospital settings



RIGID STYLET

For use with endotracheal tubes 6,0 mm and larger.

Recommended to facilitate successful intubations.



EMS BAG

Color: BLUE/RED

Size: 3.5inch/8inch



Full Cover Blades covid-19

Mac1/2/3/4/5

Option: blades with channel (Mac 3/4)





Silicone Cover Blades

Mac 1/2/3/4/5


Option: blades with channel (Mac 3/ 4)



Half Cover Blades

Mac 1/2/3/4/5

Option: blades with channel (Mac 3/ 4)



Reusable Blades

Model No.	Size(L*W*H)		Applicable patients
Miller 00		72X14.4X29.5	premature
Miller 0		84X15X29.5	newborn
Miller 1		110X15X29.5	infant
Mac 1	MacD 1	101X18X34 96X18X34	infant
Mac 2	MacD 2	115X21X36.5 115X22X50	children
Mac 3	MacD 3	130X24X40 127X22X34	adult
Mac 4	MacD 4	136X23X43 136X22X57	difficulty (obese)
Mac 5	MacD 5	150X23X51 147X22X59	difficulty



Disposable Blades

Model No.	Size(L*W*H)	Applicable patients
Mac 1	95X28X12.5	infant
Mac 2	110X29X13	children
Mac 3	120X27X13	adult
Mac 4	135X22X14	difficulty (obese)
Mac 5	120X29X16	difficulty





Equibiomedic

Venta, mantenimiento y alquiler de equipos biomédicos.



PBX: (057) (4) **444 08 38** / Cel: **310 460 7889**



www.equibiomedic.com



Medellín - Colombia

